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

## Biospin Omni Genomic DNA Extraction Kit

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
## Storage & shipping conditions

1. Store Protease K and RNase A at 2~8°C, store all other components at room temperature.
2. All components are stable for 18 months from date of receipt under proper storage condition.

## Introduction

The kit provides a simple, fast and cost effective method that can isolate high quality DNA. Using one simple protocol, high yield purified DNA can be isolated from various samples including **animal tissue**, whole blood, buffy coat, leukocytes, **bacteria**, yeast, blood spots, swabs and cultured cells. The kit utilizes silica-based membrane techniques within a convenient spin column, no need of expensive resins, toxic phenol-chloroform extractions, or time-consuming alcohol precipitation. The whole procedure takes less than 20 minutes including lysis steps, and more than 20 kb DNA could be purified by this kit. Isolated DNA can be directly applied to PCR, Southern blotting and other enzymatic reactions.

## Typical genomic DNA yields from various samples

Source	Quantity	Yield (μg)
Mammalian tissue	25mg	5-30
Mammalian blood	200μL	2-10
Bacterium	≤ 10 <sup>9</sup>	≤ 50
Cell	≤ 2×10 <sup>6</sup>	5-30
Blood spots	2-5 pieces	0.1-1

## Additional apparatus and materials required but not supplied

- |                                       |  |
|---------------------------------------|--|
| * Sterile 1.5mL microcentrifuge tubes | * 10μL/200μL/1000μL tips   |
| * Centrifuge capable of 12,000g       | * Absolute ethanol   |
| * Vortex mixer                        | * Warm bath  |
| * PBS buffer                          | * Lysozyme   |
| * Lyticase                            | * Sorbitol buffer (18.217g sorbitol + 3.7224g EDTA-2Na, dissolve into 100 mL ddH <sub>2</sub> O) |

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## Important Notes

Please add absolute ethanol to **PW Buffer** and **Wash Buffer**, mix thoroughly before use.

## Protocol

### ➤ Sample lysis

#### ■ Cell and blood

##### 1. Sample pre-processing

A. **Culture cells:** adherent cells shall be digested with trypsin, no more than  $10^7$  cells. Spin for 5min in 1500rpm, discard the supernate. Add 200μL PBS buffer into the tube, re-suspend cells.

B. **Leukocytes:** if blood volume is more than 200μL, use Red Blood Cell Lysis Buffer (Cat. # BSA06M1 or Cat. # BSA07M1) first to obtain white blood cells, add 200μL PBS buffer into the tube, re-suspend cells.

C. **Whole blood:** Add 200μl well mixed blood into tube (if the blood is less than 200μL, add to 200μL with PBS buffer).

2. Add 20μL PK solution, mix thoroughly for 15s. Add 4μL RNase A. Incubate at room temperature for 2min.
3. Add 200μL Lysis B buffer, mix thoroughly. Incubate at 56°C for 10min.
4. Add 200μL absolute ethanol, mix thoroughly.

#### ■ Animal tissue

1. Add 180μL TES Buffer and 20μL PK solution into tube.

2. Grind tissue to powder (no more than 25mg) with liquid nitrogen. Add it to the tube and mix thoroughly.

**Note:** for spleen, no more than 10mg, samples must be submerged within solution.

3. Incubate for 1-4 hours at 56°C. For samples hard to lysis, prolong lysis time or even overnight.

**Note:** vortex per 10s will facilitate lysis.

4. Centrifuge 3min at 12,000g. Carefully transfer the supernatant into a new tube.
5. Add 4μL RNase A, incubate at room temperature for 2min.
6. Add 200μL Lysis B buffer, mix thoroughly.
7. Add 200μL absolute ethanol, mix thoroughly.

#### ■ Blood spots

1. Take 2-5 pieces spots (3×3mm) into a tube.
2. Add 280μL TES Buffer and 20μL of PK solution into the tube. Mix intensively for 10s. Put it in 56°C thermostat oscillator metal bath, vortex for 60min at 900rpm.
3. Centrifuge for 3min at 12,000g. Carefully transfer the supernatant into a new tube.
4. Add 4μL RNase A, incubate at room temperature for 2min.
5. Add 300μL Lysis B buffer, mix thoroughly.
6. Add 150μL absolute ethanol, mix thoroughly.

#### ■ Gram negative bacteria lysis

1. Harvest bacterium (maximum  $2 \times 10^9$  cells) in a microcentrifuge tube by centrifuging for 1 min at maximum speed. Discard supernatant as far as possible.
2. Resuspend pellet in 180 μL TES Buffer and 20μL PK solution. Mix by pulse-vortexing intensively for 20seconds.
3. Incubate at 56°C for 30 minutes-4hours.
4. Optional: add 4μL RNase A (20mg/ml) and mix thoroughly.

5. Add 200μL LysisB buffer, and mix thoroughly.
6. Add 200μL absolute ethanol, mix thoroughly.
- **Gram positive bacteria lysis**
1. Please add Lysozyme to TET Buffer to a final concentration of 20mg / mL, mix thoroughly for 30seconds or so until the solution is clear, stored at 2-8°C.
2. Harvest bacterium (maximum 2 x 10<sup>9</sup> cells) in a microcentrifuge tube by centrifuging for 1 min at maximum speed. Discard supernatant as far as possible.
3. Resuspend pellet in 180 μL TES Buffer. Mix by pulse-vortexing intensively for 20seconds.
4. Incubate at 37°C for 30 minutes.
5. Add 20μL PK solution and mix thoroughly.
6. Add 200μL LysisB buffer, and mix thoroughly. Incubate at 55°C for 30 minutes
7. Add 200μL absolute ethanol, mix thoroughly.
- **Yeast culture lysis**
1. Resuspend the yeast cells pellet in 500 μL sorbitol buffer. Add 100~200 units lyticase and incubate at 37 °C for 30minutes. Centrifuge for 1 min at 6,000 x g. Discard supernatant.
2. Add 180μL TES Buffer and 20μL PK solution, vortex the tube for 20sec.
3. Incubate at 56°C for 30 minutes.
4. Optional: add 4μL RNase A (20mg/ml) and mix thoroughly. Incubate at room temperature for 2 minutes.
5. Add 200μL LysisB buffer, and mix thoroughly.
6. Add 200μL absolute ethanol, mix thoroughly.
- **Swab**
1. Sample preparation: Use cotton swabs in cheek wipes 10-15 times.
2. Turn the swab to a 2 mL centrifuge tube, Cut out stem part from the swab with scissors.
3. Add 500μL LysisB Buffer and 20μL PK Solution to the 2.0mL microcentrifuge tube. Mix thoroughly for 10 seconds.
4. Incubate at 56°C for 15 minutes.
5. Take 350μL the lysis product to the 1.5mL microcentrifuge tube.
6. Add 350μL of ethanol and mix thoroughly for 10 seconds.

### ➤ DNA Purification

1. Transfer all the mixture(≤750μL) to spin column. Centrifuge the mixture at 12,000×g for 1 minute. Discard flow-through.
2. Add 500μL PW buffer to the Spin column. Centrifuge the spin column at 12,000×g for 1 minute. Discard flow-through.
3. Add 500μL Wash buffer to the spin column. Centrifuge at 12,000 ×g for 1 minute. Discard flow-through. Take the spin column back to the tube.
4. Repeat step 3
5. Centrifuge the spin column at 12,000 ×g for 2 minute.
6. Place the spin column in a new 1.5 or 2.0mL micro centrifuge tube. Add 30-100μL of the Elution Buffer. Incubate at room temperature for 2 minutes. Centrifuge the mixture at 10,000×g for 1 minute. The DNA in the collection tube is ready for further analysis. If the isolated DNA sample is not going to be tested on the same day, freeze it at -20°C.

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